

Sicher – Sauber – Aufgeräumt Safe – clean – tidy

DIA-Safestore schützt DIATEST Handmessmittel und bewahrt diese sicher auf.
DIA-Safestore protects DIATEST hand measuring tools and stores them safely.





Schale für Messuhr und Anzeigeeinheit
Shell for indicator and display unit

Messmittel müssen sauber und sicher gelagert werden, damit sie ihre Messaufgabe einwandfrei und präzise erledigen können.

- Schützt DIATEST Handmessmittel und bewahrt diese sicher auf.
- Kann an jedem Messarbeitsplatz eingesetzt werden.
- Das Messgerät wird komplett montiert abgelegt, kann leicht entnommen und dadurch sofort eingesetzt werden.
- Perfekte Integration mit dem DIATEST-Messarbeitsplatz DIA-Gauge Station

For the perfect and precise completion of measuring tasks measuring instruments have to be stored in a clean and dry condition.

- Protects DIATEST hand measuring tools and stores them safely.
- Can be used at every gauge station.
- The measuring instrument is stored completely assembled, can easily be removed and is thus immediately ready for use.
- Perfect integration with the DIA-Gauge Station

Perfekte Aufbewahrung von:

- DIATRON1000
- Bohrungsmessdorne (Ø 7,0 – 40,0 mm)
- IKT, AKT (Außen-Ø max. 50,0 mm), KT, KT-B
- TD-Gauges+ (DI)
- Gauge-2Dimensions
- Tastköpfe
- Messköpfe (kleine Größen)
- DIN 2250-C Einstellringe (Ø 5,0 – 40,0 mm)

Perfect storage of:

- DIATRON1000
- BMD Plug Gauges (Ø 7 – 40 mm)
- IKT, AKT (O.D.-Ø max. 50 mm), KT, KT-B
- TD-Gauges+ (DI)
- Gauge-2Dimensions
- Split-Ball Probes
- Plunger Probes (small sizes)
- DIN 2250-C Setting rings (Ø 5,0 – 40,0 mm)



Kombination verschiedener Ablagen
Combination of different storages



Ablage für TD-Gauge+ (DI)
Storage for TD-Gauge+ (DI)

➤ Für ein Angebot wenden Sie sich bitte an Ihren DIATEST-Händler vor Ort. ➤ Please contact your local DIATEST dealer for an offer.

DIATEST Hermann Költgen GmbH

Schottener Weg 6
D-64289 Darmstadt
Germany

☎ +49 61519790
☎ +49 6151979111
✉ info@diatest.com
www.diatest.com

